

Saliva Tracer 7

Australia
Breathalysers

✓ Verified to AS/NZS 4760:2019

✓ 7 x Different Drug Groups

✓ Latest Technology

Ideal Application: Accelerated Outcomes Through The Use of Advanced Technology

The **DrugCheckExpress® SalivaTracer™ 7** uses advanced saliva testing technology to deliver the fastest screening speed available. It provides negative results in under 2 minutes, offering a simple, hygienic process that supports rapid decision-making and enables streamlined testing procedures.



CORE FEATURES

Advanced Saliva Testing Technology:

Built-In Saturation Indicator turns red once sufficient saliva has been collected to support accurate test activation and reduce collection errors.

Rapid Turnaround Results:

Negative results can be read in as little as 2 minutes for fast and efficient on-site screening.

Precision Collection Swab:

Saliva-based testing offers a convenient alternative to urine collection with reduced privacy concerns during sample collection.

Modern Low-Profile Design:

Simple collection and interpretation process supports efficient workplace drug screening procedures.

High-Efficiency Screening Solution:

Saliva testing is designed to help identify recent drug use and potential impairment at the time of testing.

TECHNICAL SPECS

- **Drugs Detected:** Cannabis/THC (THC) Methamphetamines (MET) Cocaine (COC) Amphetamines (AMP) Benzodiazepines (BZO) Opiates (OPI) Oxycodone (OXY)
- **Cut-Off Levels:** Cannabis/THC (THC): 15 ng/mL
Methamphetamines (MET): 50 ng/mL Cocaine (COC): 50 ng/mL
Amphetamines (AMP): 50 ng/mL Opiates (OPI): 50 ng/mL
Oxycodone (OXY): 40 ng/mL Benzodiazepines (BZO): 10 ng/mL
- **Units per Box:** 25
- **Sample Type:** Saliva
- **Dimensions:** 13.9 x 3.7 x 0.7 cm
- **Weight:** 18g

