

✓ Verified to AS/NZS 4308:2023

✓ 6 x Different Drug Groups

✓ Leak Resistant Key Security

Ideal Application: Longer Detection Window with Key Security

The **ToxTest** KC 600 Split Key Urine Drug Test Cup combines advanced testing technology with a secure split key design for reliable on-site urine screening. Featuring a longer detection window, secure leak-resistant design, and easy sample handling, it supports accurate and efficient workplace testing.



CORE FEATURES

Secure Split-Cup System:

Integrated split specimen design allows for additional laboratory confirmation testing if required.

Tamper-Resistant Key Security:

Includes checks for pH, nitrite, oxidants, specific gravity, and creatinine to help identify compromised samples.

Extended Detection Capability:

Urine-based testing supports a longer detection period compared to saliva testing, making it suitable for workplace, random, and pre-employment screening.

Built-In Temperature Monitoring:

It has an inbuilt temperature strip that lets you monitor the temperature of the urine sample so you can ensure that the sample hasn't been tampered with.

Secure Sample Integrity:

Simple cup-based format with clear result windows supports convenient sample collection and interpretation.

TECHNICAL SPECS

- **Drugs Detected:** Cannabis/THC (THC) Methamphetamines (MET) Cocaine (COC) Amphetamines (AMP) Benzodiazepines (BZO) Opiates (OPI)
- **Cut-Off Levels:** Cannabis/THC (THC): 50 ng/mL
Methamphetamines (MET): 300 ng/mL Cocaine (COC): 150 ng/mL
Amphetamines (AMP): 300 ng/mL Opiates (OPI): 300 ng/mL
Benzodiazepines (BZO): 200 ng/mL
- **Units per Box:** 25
- **Sample Type:** Urine
- **Adulterants Detected:** Creatinine - pH
Specific Gravity - Oxidants - Nitrite
- **Temperature Strip:** Yes
- **Dimensions:** 8.7 x 6.4 x 7.4 cm
- **Weight:** 69g

